


Try It!  
You'll Like It!

# ORTHOMERICA'S NEW FOAM HUMERAL FRACTURE BRACE

*Compression Molded Foam*

- ☆ **Ultralight compression molded foam provides total contact support**
- ☆ **Improved strapping system eases application and removal process**
- ☆ **Dual tongue system maximizes compression and fit**
- ☆ **Comfortable interface with skin for increased patient compliance**
- ☆ **"Shoulder Cap" limits migration yet allows full shoulder movement**
- ☆ **Removable distal strap allows orthotist to optimize positioning to accommodate length**
- ☆ **Easily trimmed with scissors for clinical applications**
- ☆ **Universal design in five sizes fits right or left arm**



The Compression molded foam humeral fracture brace effectively stabilizes diaphyseal fractures to provide stability and reduce pain. Lightweight, strong, and comfortable, this orthosis provides an intimate fit for effective fracture management and compliance.

SIZE	MODEL NO.	DISTAL CIRC.	PROXIMAL CIRC.	BRACE LENGTH
X-Small	1229	6"-8"	8"-10"	11"
Small	1230	7"-9"	9"-11"	12"
Medium	1231	9"-11"	11"-13"	13"
Large	1232	11"-13"	13"-15"	13"
X-Large	1233	13"-15"	15"-17"	14"

*Kit includes 2 socks*

*Unique dual tongues molded into the medial component can be trimmed to provide compression*

