



STAR Family of Cranial Remolding Orthoses

General Information

Used to Treat*

Measurements

Designer

Manufacturer

Materials

Style



**STARband™
Side-Opening Band**

Most widely prescribed headband design

Plagiocephaly
Brachycephaly
Scaphocephaly
Post-operative Cranio-synostosis

Custom fabricated from a plaster cast or STAR scan

Orthomerica

Exclusively fabricated by Orthomerica

Plastic shell with foam liner

Side and top openings



**STARband™
Bi-Valve**

Most used for scaphocephaly head shapes

Plagiocephaly
Brachycephaly
Scaphocephaly
Post-operative Cranio-synostosis

Custom fabricated from a plaster cast or STAR scan

Orthomerica

Exclusively fabricated by Orthomerica

Plastic shells with foam liner

Overlap design



**STARlight™
Side Opening Band**

Clear material allows visual inspection

Plagiocephaly
Brachycephaly
Scaphocephaly
Post-operative Cranio-synostosis

Custom fabricated from a plaster cast or STAR scan

Orthomerica

Exclusively fabricated by Orthomerica

Clear plastic shell with foam pads

Side and top openings



**STARlight™
Bi-valve**

Clear material allows visual inspection

Plagiocephaly
Brachycephaly
Scaphocephaly
Post-operative Cranio-synostosis

Custom fabricated from a plaster cast or STAR scan

Orthomerica

Exclusively fabricated by Orthomerica

Clear plastic shells with foam pads

Overlap design



**STARlight™ PRO
Side Opening Band**

First helmet for post-endoscopic craniectomy

Post-operative Cranio-synostosis release

Custom fabricated from a plaster cast or STAR scan

Dr. Jimenez & Dr. Barone

Exclusively fabricated by Orthomerica

Clear plastic shell with foam pads

One piece living hinge



**Clarren
Helmet**

First helmet used to treat head shape deformities

Plagiocephaly
Brachycephaly
Scaphocephaly

Custom fabricated from a plaster cast or STAR scan

Dr. Sterling Clarren

Exclusively fabricated by Orthomerica


Plastic shell with foam liner

Full plastic helmet with chin strap

*Selection of cranial remolding orthosis design and treatment protocol is subject to clinical evaluation and physician prescription.



Before Treatment



After STARband Treatment

ORTHOMERICA®
877-737-8444 | orthomerica.com