



UPPER EXTREMITY CUSTOM-TO-CAST



MEASUREMENT LEGEND

- A = Circumference at Proximal Upper Arm
- B = Circumference at Distal Upper Arm
- C = Circumference at Proximal Forearm
- D = Circumference at Distal Forearm
- E₁ = Length from Axilla to Elbow Center
- E₂ = Length from Acromion Process to Elbow Center
- F = Length from Elbow Center to Mid Palm
- G = Length from Elbow Center to Styloid
- H = Caliper M/L at Elbow
- J = Palm Width



Extend Over Shoulder

Yes No

Plastic: _____

Liner: _____

UPGRADES

- Upgrade to Bilateral Adjustable Elbow Joints
(Standard includes adjustable lateral, non-adjustable medial joint.)
- Upgrade to Bilateral R.O.M. Wrist Joint
(adjustable lateral and medial joints.)
- Other _____

ORDER INFORMATION

Date: _____ P.O.#: _____
 Facility to be Billed: _____
 Ship to Address: _____

 Phone: _____ Fax: _____
 Date Required: _____
 Ship Via: _____ On (Date) _____
 Contact: _____

PATIENT INFORMATION

Date: _____
 Patient Name: _____
 Age: _____ Sex: _____ Height: _____ Weight: _____
 Diagnosis: _____

Includes 2 socks and a sling

**Fax This Form to 877.737.8445
 or Call 877.737.8444**