

# Newport<sup>®</sup> lite<sup>™</sup>

The new Newport lite hip orthosis from Orthomerica is designed to provide protection and support for patients at risk for problems related to hip instability.

## INDICATIONS

- Arthroscopic hip procedures
- Labral repair, with or without gluteus medius repair
- Post-op proximal hamstring repair
- Patients at risk to dislocate

## CONTRAINDICATIONS

- Global instability
- Patients at risk for an anterior dislocation
- Chronic dislocation

## FEATURES/BENEFITS

- Protects the hip joint capsule by providing a limited and safe range of motion following surgical procedures
- Well-tolerated by patients
- Virtual joint is easily set to control hip flexion, extension, adduction and abduction
- Proprioceptive reminder
- Universal for right or left side
- Easily customized to provide an exact fit
- One size fits most
- Laterally placed abduction paddle helps to maintain the head of the femur in the acetabulum
- Compression belt design with pull tabs anchors the hip joint and thigh cuff to promote optimal alignment
- Accommodates hip development through ball joint mechanism





#### Compression Belt

Dual pull-tab system is simple for the patient to tighten as well as don and doff.

**Three Point Design**  
Helps push the femoral head into the acetabulum.



Newport® lite™ Hip Orthosis *one size fits most waists 26"–66" (66cm–167cm)*

Item	Model #
Complete Orthosis	3662
Abduction Control Joint	3662.03
Pelvic Component	3662.04
Replacement Pad Kit	3662.05

6333 North Orange Blossom Trail • Orlando FL 32810 • USA

**ORTHOMERICA®**

800-446-6770 | [www.orthomerica.com](http://www.orthomerica.com)



Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands