



Newport[®] 3

HIP SYSTEM

The Newport 3 Hip Orthosis from Orthomerica has been the industry standard for over a decade. Innovations like the patented clamping disks, allowing precise contouring of the pelvic section and the Universal Virtual Plus Joint, which provides unprecedented strength in a small component, have made the Newport Hip Orthosis the most popular choice for physicians, orthotists, and patients.

Features

- Pelvic section is designed for either right or left hip joint attachment without disassembling the orthosis
- Patented locking disks adjust for hip development and circumference
- Hip joint can be used on right or left side and is adjustable for flexion, extension, abduction, and adduction
- Low profile thigh component clears the groin and does not impinge on knee flexion
- Newport hip orthosis has been proven to prevent dislocation in a scientific study (JPO, Vol. 6, Number 1, Winter 1994.)

Indications

- Post-operative hip revision patients
- Primary arthroplasty patients at risk to dislocate
- Patients needing stability after dislocation
- Inoperable patients requiring hip stabilization
- Patients who can benefit from a hip orthosis to reinforce hip precautions

* *Standard Newport 3 Hip System*

* *Freight carrier to be specified by Orthomerica. Freight special only applies to prefabricated Newport 3 & 4 Hip Systems shipped out of the Orlando facility.*