

Orthomerica's STAR Family of Cranial Remolding Orthoses

General Information

Used to Treat*

Measurements

Designer

Materials

Style



STARband™ Side-Opening Band



STARband™ Bi-Valve



STARlight™ Side Opening Band



STARlight™ Bi-valve



STARlight™ PRO Side Opening Band

Most widely prescribed headband design

Most used for scaphocephaly head shapes

Clear material allows visual inspection

Clear material allows visual inspection

For post-endoscopic craniectomy

Plagiocephaly
Brachycephaly
Scaphocephaly
Post-operative Cranio-synostosis

Plagiocephaly
Brachycephaly
Scaphocephaly
Post-operative Cranio-synostosis

Plagiocephaly
Brachycephaly
Scaphocephaly
Post-operative Cranio-synostosis

Plagiocephaly
Brachycephaly
Scaphocephaly
Post-operative Cranio-synostosis

Post-operative Cranio-synostosis release

Custom fabricated from a plaster cast or scan

Custom fabricated from a plaster cast or scan

Custom fabricated from a plaster cast or scan

Custom fabricated from a plaster cast or scan

Custom fabricated from a plaster cast or scan

Orthomerica

Orthomerica

Orthomerica

Orthomerica

Dr. Jimenez & Dr. Barone

Plastic shell with foam liner

Plastic shells with foam liner

Clear plastic shell with foam pads

Clear plastic shells with foam pads

Clear plastic shell with foam pads

Side and top openings

Overlap design

Side and top openings

Overlap design

One piece living hinge

*Selection of cranial remolding orthosis design and treatment protocol is subject to clinical evaluation and physician prescription.



Before Treatment



After STARband Treatment

ORTHOMERICA[®]
6333 North Orange Blossom Trail • Orlando FL 32810 • USA
877-737-8444 | www.orthomerica.com

EC REP Emergo Europe
Molenstraat 15, 2513 BH,
The Hague, The Netherlands