

Orthomerica's STAR Family of Cranial Remolding Orthoses



General Information

Used to Treat*

Measurements

Designer

Materials

Style



STARband™ Side-Opening Band




STARband™ Bi-Valve



STARlight™ Side Opening Band



STARlight™ Bi-valve



STARlight™ PRO Side Opening Band

Most widely prescribed headband design

Plagiocephaly
Brachycephaly
Scaphocephaly
Post-operative Cranio-synostosis

Custom fabricated from a plaster cast or scan

Orthomerica

Plastic shell with foam liner

Side and top openings

Most used for scaphocephaly head shapes

Plagiocephaly
Brachycephaly
Scaphocephaly
Post-operative Cranio-synostosis

Custom fabricated from a plaster cast or scan

Orthomerica

Plastic shells with foam liner

Overlap design

Clear material allows visual inspection

Plagiocephaly
Brachycephaly
Scaphocephaly
Post-operative Cranio-synostosis

Custom fabricated from a plaster cast or scan

Orthomerica

Clear plastic shell with foam pads

Side and top openings

Clear material allows visual inspection

Plagiocephaly
Brachycephaly
Scaphocephaly
Post-operative Cranio-synostosis

Custom fabricated from a plaster cast or scan

Orthomerica

Clear plastic shells with foam pads

Overlap design

For post-endoscopic craniectomy

Post-operative Cranio-synostosis release

Custom fabricated from a plaster cast or scan

Dr. Jimenez & Dr. Barone

Clear plastic shell with foam pads

One piece living hinge

*Selection of cranial remolding orthosis design and treatment protocol is subject to clinical evaluation and physician prescription.



Before Treatment



After STARband Treatment

ORTHOMERICA[®]
6333 North Orange Blossom Trail • Orlando FL 32810 • USA
877-737-8444 | www.orthomerica.com

EC REP Emergo Europe
Molenstraat 15, 2513 BH,
The Hague, The Netherlands